

***TherAdapt<sup>®</sup> Products Inc.***  
(800) 261-4919 phone (866) 892-2478 fax

**Adult Angle/Straddle  
Bench**

**www.theradapt.com**

**Please read all of the following information before using this product.**

*TherAdapt<sup>®</sup> Products Inc. has designed and manufactured this product with the well being of the user and caregiver(s) in mind. It is our endeavor to provide products which are attractive, comfortable, useful, highly adjustable, long lasting and safe. As when using any product, precautions must be taken.*



**WARNING:** This product must only be used under the recommendation of a physician and / or licensed therapist and with adult supervision. Never leave user unattended. It is your responsibility to properly use and maintain your Adult Angle Straddle Bench. Failure to follow USE and CARE instructions or ASSEMBLE the bench properly as outlined may result in bodily injury.

**THE MANUFACTURER IS NOT RESPONSIBLE FOR INJURIES RESULTING FROM ANY CUSTOMIZATION OR MODIFICATION OF THE EQUIPMENT THAT WAS NOT COMPLETED BY THE MANUFACTURER.**

**USE:** The Adult Angle Straddle Bench is designed for a single purpose: as a freestanding bench. It is to be used by a single individual that is seated in the center of the bench with their feet firmly placed on the floor. **No other use is permitted.**

**CARE:** Take care to inspect the ongoing condition of all parts of your Adult Angle Straddle Bench. Verify that all knobs, bolts, and washers are properly positioned and secure. Contact the manufacturer for replacement or repair of necessary items. All components may be cleaned with mild soap and water, or a gentle sanitizing spray. **Do not** use abrasive cleaners. Wipe wood components dry after cleaning. Excess amounts of liquid on the wood may damage the finish. Always store indoors.

**WARRANTY:** The manufacturer provides a full one (1) year warranty, beginning the date shipped by the manufacturer. The warranty provides full replacement of all parts against breakage when used in compliance with the manufacturer's USE and CARE instructions. Use or modification beyond the normal capacity or range of adaptability designed into the products will void the warranty. The foregoing is TherAdapt's sole liability for consequential damages.

**PLEASE COMPLETE AND RETAIN FOR YOUR RECORDS**

**Adult Angle Straddle Bench AASB-\_\_\_\_\_ Serial Number: \_\_\_\_\_**

**Purchase date: \_\_\_\_\_**

## *TherAdapt*<sup>®</sup> Adult Angle/Straddle Bench (AASB-500) Instructions

### Packing List:

- (1) Bench top component with suede Naugahyde<sup>®</sup> padding
- (2) Rectangular side pieces with (4) ¼ -20 counter-sunk bolts each
- (2) Curved rectangular floor pieces
- (2) Curved rocker bottom pieces
- (8) Thru-knobs with ¼ -20 threaded inserts
- (8) Nylon washers



### 1. Attach the curved, rectangular floor pieces.

Remove the (4) lower thru-knobs, attach floor pieces, and replace thru-knobs. **Tighten knobs securely.**

### 2. Adjust the height.

To adjust the height, we recommend standing the bench on one end, and removing the pair of knobs on the inner surface of the upper rectangular sidepiece. The combined height and angle adjustable pieces can be lifted off and the appropriate height determined by putting bolts through one of the four pairs of holes drilled in the sides.



### 3. Adjust the angle.

The outermost, curved, rectangular sidepieces determine the angle setting. Remove the outer pair of knobs and lift off the angle piece. **Note:** This piece must be lifted off completely to make adjustments. The bench angles from a horizontal setting using the bottom hole to a 10°, 15°, 20°, 25°, or 30° angle using the five ascending holes.

### 4. Convert the angle bench to a straddle bench.

Remove the curved pieces from their storage position under the bench top. Remove the curved-rectangular pieces and replace with the rocker bottoms. Adjust the height as above.

**Tighten all knobs securely before using.**



**FOR PROFESSIONAL ASSISTANCE CALL 1-800-261-4919**