



Adaptive School Chairs are designed to provide optimum seating posture for participation in classroom, while eating and during therapy.







The First Class School Chair is ideal for children with mild to more involved physical challenges for use in classrooms, clinics or homes. Its adjustablility and wide range of accessories meet the needs of each child as they grow and develop.

HCPCS E1399

**Tool Free Adjustment!** Simplifies use for multiple children in a classroom and clinic setting.

### First Class School Chair

Item #	Product Size	UOM
FC 2000N	Small	1/bx
FC 4000N	Large	1/bx

- Seat height and depth adjustable solid seat with a 15° anterior or posterior tilt.
- Armrests provide lateral limitations as well as a stable resting surface for arms. Flip-up armrests (Figure A) allow easy transfers and they are height adjustable. The tray accessory mounts onto the armrest.
- Padded Hip Belt 2-point hip belt provides pelvic support and maintains proper hip positioning.

anacifications	E0 0000N	FO 4000N
specifications	FC 2000N	FC 4000N
Seat Width	12"	16"
Seat Depth	9.5" - 14.5"	13" - 17"
Seat Height to Floor	10" - 14"	15" - 23"
Seat Height to Footplate	9" - 12.5"	9.5" - 17.5"
Back Height	15"	17"
Length of Armrest	10"	10"
Seat-to-Armrest Height	5" - 10"	4" - 13.5"
Distance Between Armrests	13.5"	17.5"
Hip Guides	3" H x 10" D	3" H x 10" D



### **Multi-Axis Headrest**





**HCPCS** 

E1399

Item #	For Use With	UOM
FC 8000N	FC2000N & FC4000N	1/bx

- · Height, depth & angle adjustable.
- Ball and socket style mounting allows for unlimited pad positioning.
- · Long vertical rod provides extensive height adjustment capability.
- Contoured pad rotates and can be placed off-center.
- Adaptable to any solid back wheelchair.

# **Mobility Legs**





**HCPCS** E1399

Item #	For Use With	UOM
FC 2029N	FC 2000N	1set/bx
FC 4029N	FC 4000N	1set/bx

- Height adjustable legs with wheels offer safe mobility.
- · Rear legs include brakes and anti-tippers.

#### **Footrests**





Item #	For Use With	UOM	
FC 2028	FC 2000N	1pr/bx	
FC 4028	FC 4000N	1pr/bx	

• Two-piece flip-up footplate with heel loop.

### **Leg Extensions**



Item #	For Use With	UOM
FC 8050	FC 2000N & FC 4000N	Set of 4/bx

• Set of four leg extensions for custom height.



### **Anti-tippers**



1Page
HCPC <b>E139</b>

Item #	For Use With	UOM
FC 2030N	FC 2000N	1/bx
FC 4030N	FC 4000N	1/bx

• Anti-tippers provide additional stability and prevent chair from tipping backwards.

# **Hip Guides**





E1399

Item #	For Use With	UOM
FC 8027N	FC 2000N & FC 4000N	1pr/bx

- · Width and depth adjustable.
- Supports user's pelvis and prevents adduction.

### **Tray**

FC 4024



1PACK
HCPCS
E1399

1/bx

Item #	For Use With	UOM
FC 2024	FC 2000N	1/bx

- FC 4000N • Tray attaches onto armrest and is height adjustable.
- · Provides trunk stabilization and an easy-to-clean surface with a safety lip.
- · Ideal for school, play and feeding.

### **Push Handle**



Item #	For Use With	UOM
FC 8001	FC 2000N & FC 4000N	1pr/bx

· For use with the First Class chair with the Mobility legs, the push handles allows for easy maneuverability and transportation.











## Support Kit

Item #	Product Size	UOM
FC 8025N	Small	1set/bx
EC 8025 I N	Large	1set/hy

- Support Kit includes Trunk Harness, Abductor and Lateral Supports.
- Trunk Harness:
- Adjustable zippered trunk harness provides focused support on the sternum and rib cage.
- its semi-stretch material allows controlled movement, prevents forward slumping and also maintains correct midline seated position.
- Straps with cam lock buckles secure the trunk harness.
- Abductor:
- · Depth adjustable
- Positions user's knees controls spastic muscles that might cause their legs to cross.
- · Lateral Supports:
- Assists proper stabilization and alignment of the upper body.
- Height and width adjustable.
- Sold as a pair and each lateral is individually adjustable.

### **Trunk Harness**

Included with Support Kit, can be purchased on its own

Item #	<b>Product Size</b>	UOM
VTH 001N	Trunk Harness - Small	1/bx
VTH 001 LN	Trunk Harness - Large	1/bx

- Adjustable zippered trunk harness provides focused support on the sternum and rib cage.
- Its semi-stretch material allows controlled movement, prevents forward slumping and also maintains correct midline seated position.
- Straps with cam lock buckles secure the trunk harness.

specifications	VTH 001	VTH 001L
1A	8"	13"
2A	12"	14"
3A	10.5"	13"
Strap Length	12"	16"



