

MSS Tilt & Recline

The MSS Tilt & Recline is designed for children with moderate to severe physical involvement. It features an easy tilt and recline mechanism with up to 45 degrees tilt and the back reclines from 80° up to 140°. The basic model includes all the features listed below.

MS 2000 Small, 1/cs

MS 3000 Medium, 1/cs

FEATURES:

- Footrest – adjustable in height and depth.
- Seat – adjustable in depth.
- Armrest – adjustable in height and angle.
- Back – adjustable in height and angle.
- Headrest – adjustable in height, depth and angle. It rotates and can be placed off center.
- Lateral supports – individually adjustable in height and width.
- Abductor – adjustable in height and depth. It can be placed between the thighs or turned around and placed between the knees.
- Belt – adjustable pelvic belt.
- Tray – horizontal positioning device.
- Harness – zippered trunk harness in an H configuration for anterior chest support.
- Cushion covers are removable and washable.
- Two year limited warranty.

BASE:

MS 0014 Hi Base for MSS Tilt & Recline, 1/cs

MS 0015 Low Base for MSS Tilt & Recline, 1/cs

Wheeled mobility bases fit both sizes MSS Tilt & Recline chairs. It allows the MSS Tilt & Recline to be used as floor sitter or a high chair. The Hi Base will bring the seating system up to table height.



DIMENSIONS	MS 2000	MS3000
SEAT WIDTH	9"	10.5"
SEAT DEPTH	5" - 14"	5" - 14"
BACK HEIGHT	10" - 12.5"	12.5" - 16"
SEAT TO FOOTREST HEIGHT	5" - 10.5"	5" - 10.5"
FOOTREST DEPTH	7.5"	7.5"
ARMREST TO SEAT HEIGHT	5" - 9"	5" - 9"
FLOOR TO SEAT HEIGHT		
AS FLOOR SITTER	13.5"	13.5"
WITH LOW BASE	16"	16"
WITH HIGH BASE	22.5"	22.5"
SPACE BETWEEN LATERALS	7" - 11"	7" - 11"
HEADREST HEIGHT ADJUSTMENT	Adds 3" - 7"	Adds 3" - 7"
WEIGHT LIMIT	75 lbs.	100 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

positioning
seating

Wenzel[®] Re/hab

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM



444

MSS TILT & RECLINE ACCESSORIES:

Headrest Wedges

WD 8000 Wedges, 1pr/bx

Wedges can be placed anywhere on the headrest for additional positioning.



Vinyl Accessories

MS 2018 Complete Set of Vinyl Covers for MS 2000, 1pr/bx

MS 3018 Complete Set of Vinyl Covers for MS 3000, 1pr/bx

Clear Tray

MS 0023 Tray for MS 2000 & MS 3000, 1/bx

Clear Horizontal positioning device.

Booster Kit

MS 4000 Booster Kit for MSS Tilt & Recline insert, 1/bx

The booster kit, consisting of a pad and straps, can be attached to the MSS insert and easily converts it to a booster chair. It can then be secured to most standard chairs and used at dining room and kitchen tables.



Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

MSS TILT & RECLINE ACCESSORIES:

Traveler Stroller Base

MS 5000 Traveler Base, 1/cs

A lightweight stroller base accommodates the MSS Tilt & Recline seating system. The MSS Tilt & Recline seating system is easily removed from its side frames and inserted into the Traveler frame.

FEATURES:

- Front swivel wheels.
- Oversized rear wheels.
- Retains all the Mss Tilt & Recline's positioning functions.
- Three tilt positions.
- Folds easily with the seating system in place.
- 2 year limited warranty.



TRAVELER ACCESSORIES:

Canopy

MS 5510 Canopy, 1/bx

Basket

MS 5520 Storage Basket, 1/bx



positioning
seating accessories

Wenzel *lite* Re/hab
To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.