

FCI Walker

The FCI Walker Gait Trainer is a hands-free device that encourages upright mobility and promotes independence.

The user's upright position, with a slight forward lean, encourages and stimulates forward stepping. This device can be used in two distinct ways:

Anterior - Users are forward facing for upper extremity and/or trunk support while using the attached tray as a surface for other activities or meals when stationary.

Posterior - Users face backwards to achieve closer access to involved play: to kick a ball, open doors or simply to better explore their environment.



The frame height, pommel height, chest pad height and tray heights are each individually adjustable.

The posterior support component can be removed to facilitate transference of the child in or out of the walker.

The chest support can be adjusted to change the angle of body (trunk) to a more forward prone position if either more support is needed and/or to facilitate forward movement such as stepping.

Included Features

- ▶ Clear plastic tray
- ▶ 3" swivel castors
- ▶ Back pad (S)
- ▶ Chest pad (S)
- ▶ 7" Pommel on Mini
- ▶ 9" Pommel on JR & SR

Optional Features

- ▶ 5" castors (3" Standard)
- ▶ Wheel locks (Brakes)
- ▶ Variety of directionals available
- ▶ High back customized to supplied dimensions; typically three-sided; extend to shoulders or to top of user's head
- ▶ Pelvic wedge (various sizes)
- ▶ Head rest (with hardware)
- ▶ Abductor pommel



MINI-POSTERIOR USE



MINI-ANTERIOR USE



JR



SR