# SNUG SERT











Mobile standing frame for active children



# SNUG SEAT

The new Rabbit is available in 4 sizes and suitable for children/teenagers between 1 and 18 years. The functionality of the Rabbit stays the same. It still provides the opportunity to explore and enjoy new and exciting adventures. The standing position entices the user to actively participate in things around him or her. The mobility of the Rabbit is similar to that experienced

by active wheelchair users.













Foot plate with heel stops

Knee support, hinged









ral	hh	iŀ	тм	bu	R	82°	
		Ш		by			

Measurements				
	Rabbit size 1 w/22" wheels cm (inch)	Rabbit size 2 w/26" wheels cm (inch)	Rabbit size 3 w/32" wheels cm (inch)	Rabbit size 4 w/36" wheels cm (inch)
Length, frame (L)	62 (241/4")	75 (29 <sup>1</sup> / <sub>4</sub> ")	75 (291/4")	85 (331/4")
Length, total	62 (24 <sup>1</sup> / <sub>4</sub> ")	75 (29 <sup>1</sup> / <sub>4</sub> ")	80 (311/4")	90 (35")
Width (W)	61½ (24″)	68 (26 <sup>1</sup> / <sub>2</sub> ")	69 (27")	70 (271/4")
Height from floor. Frame in highest position (H1)	100 (39")	118 (46")	131 (51")	153 (59 <sup>3</sup> / <sub>4</sub> ")
Height from floor. Frame in lowest position (H1)	67 (261/4")	85 (33 <sup>1</sup> / <sub>4</sub> ")	99 (381/2")	120 (463/4")
Height from foot plate. Frame in highest position (H2)	74 (29")	108 (42")	119 (46 <sup>1</sup> / <sub>2</sub> ")	139 (54 <sup>1</sup> / <sub>4</sub> ")
Height from foot plate. Frame in lowest position (H2)	57 (221/4")	75 (29 <sup>1</sup> / <sub>4</sub> ")	86 (331/2")	106 (411/4")
Angle of the centre bar (A1)	0 to +30°	0 to +30°	0 to +30°	0 to +30°
Angle of the foot plate (A2)	+15° to -10°	+15° to -10°	+15° to -10°	+15° to -10°
Camber	8°	4º	4º	4º
Weight	13,5 kg (29,5 lb)	15,5 kg (34 lb)	16 kg (35 lb)	18,5 kg (40,5 lb)
Max. load/user weight	40 kg (88 lb)	50 kg (110 lb)	60 kg (132 lb)	70 kg (154 lb)

# Snug Seat Inc.

12801 E. Independence Blvd. P. O. Box 1739 Matthews NC 28106 USA Tel 1 800 336 7684

Fax 704 882 0751

E-mail: information@snugseat.com Web: www.snugseat.com



# SNUG SEAT



# rabbit™

#### Angle adjustable center bar

The Rabbit has a simple, easy to use and controlled angle adjustment, making it possible to obtain the best standing position.

Central brake User operated brake system

which locks the wheels for

steady and secure transfer.



Positioning the footplate closer to the floor ensures easy placement and positioning of the child. The child

The standing position enables the child to be an active

participant, while stimulating senses and giving the

child more independence. The Rabbit is secure and re-

liable and easy to adjust, meeting the needs of the in-

dividual child. It provides children with the opportunity

to explore and enjoy new and exciting adventures. An

accessible and simple design ensures the user maxi-

mum mobility in the standing position.

# is at the same level as their play mates.













Quick release wheels

The two driving wheels are equipped with quick release for easy mounting. The wheels are available in a variety of sizes for each size Rabbit



Prepared for accessories

All sizes of Rabbit can be equipped with all sizes of accessories. The Rabbit is designed with slots in the center bar for mounting the accessories





rabbit by RBZ 6 81 80 nothing compares to a smile!





**STANDARD CONFIGURATION** 

#### Standard

As standard the Rabbit is delivered with frame and centre bar. Please order accessories from the following pages.



Frame colors
Available in different colors for size1, 2+3, 4.





for size 4 for size 2+3

Driving wheel
Please choose size of driving
wheel to complete the configuration of your Rabbit.






for size 2+3



**MEASUREMENTS** 





Measurements				
	Rabbit size 1 w/22" wheels cm (inch)	Rabbit size 2 w/26" wheels cm (inch)	Rabbit size 3 w/32" wheels cm (inch)	Rabbit size 4 w/36" wheels cm (inch)
Length, frame (L)	62 (241/4")	75 (291/4")	75 (291/4")	85 (331/4")
Length, total	62 (241/4")	75 (291/4")	80 (311/4")	90 (35")
Length, centre bar	50 (191/2")	70 (271/4")	80 (311/4")	100 (39")
Width (W)	611/2 (24")	68 (261/2")	69 (27")	70 (271/4")
Height from floor (H1) *	67-100 (261/4-39")	85-118 (331/4-46")	99-131 (381/2-51")	120-153 (463/4-593/4")
Height from foot plate (H2) *	55-74 (211/2-29")	75-108 (291/4-42")	83-119 (321/2-461/2")	96-142 (371/2-551/2")
Angle of the centre bar (A1)	0 to +30°	0 to +30°	0 to +30°	0 to +30°
Angle of the foot plate (A2)	+15° to -10°	+15° to -10°	+15° to -10°	+15° to -10°
Camber	80	6°	4º	4º
Weight	13,5 kg (29,5 lb)	15,5 kg (34 lb)	16 kg (35 lb)	18,5 kg (40,5 lb)
Max. load/user weight	40 kg (88 lb)	50 kg (110 lb)	60 kg (132 lb)	70 kg (154 lb)

\* The Rabbit is measured in upright position and with the foot plate in lowest position. For smaller children it is possible to raise the foot plate.

See Rabbit and Toucan accessories on pages 70 and 71.



SNUG SERT



# SNUG SEAT

#### ACCESSORIES

#### Foot plate Ensure the correct position of the feet.

# size measurements

W:12"xD:12" W:17"xD:14"

4+5 W:17"xD:14"

# Hinged sandals

Ensure easy and safe transfer from e.g. wheelchairs. Shown on individual footplates - also available on one-piece footplate.

#### size measurements

W:33/4"xD:8"

2+3 W:4½"xD:9½"

4+5 W:6"xD:12"

## Knee supports, hinged

Position the legs correctly while spreading the pressure above and below the knee.

#### size measurements

W:3" 2+3 W:33/a'

4+5 W:41/2"

### Abdominal support Support the child's trunk

position. The abdominal support pad can be turned upside down.

size measurements

one H:5"xW:71/2"



#### Heel stops

Ensure the correct position of the feet

size measurements

H:11/2"xW:4" H:4"xW:53/4"

## Duck feet

Ensure the correct position

size measurements

1+2 W:4" x D:81/4" 3+4 W:41/2" x D:9"

## Knee supports Position the legs correctly.

measurements

H:21/5"xW:21/5" 2+3 H:23/4"xW:31/4"

4+5 H:23/4"xW:4"

#### Hip support, front For additional support

and stability. The hip support pads can be turned upside

measurements

H:41/4"xW:6"xD:3"

2+3 H:41/4"xW:8"-101/2"xD:3"

4+5 H:61/4"xW:10"-111/2"xD:4"

#### Hip support, rear For additional support and

stability. The hip support pads can be turned upside down.

size measurements

4+5 H:61/4"xW:10"-111/2"xD:4"



#### Distance piece

To eliminate the gap between the two chest support pads.

size measurements

S H: 11 cm

H: 16 cm



### Back support

Chest support

Model Pony/Bronco.

For support of the upper body.

size measurements

1 H:14xW:21xD:12

2 H:14xW:24xD:12

3 H:14xW:29xD:12

Stabilizes the pelvis from behind in the prone or upright standing position.

size measurements

one H:13 x W:19 cm



#### Chest support

The pads can be turned upside down. Size 1 is available with either one or two pieces pads other sizes are one piece pads.

size measurements

H:11xW:15½xD:7½

H:11xW:20-27xD:71/2 2+3 H:11xW:20-27xD:71/2

4+5 H:16xW:25-32xD:10<sup>1</sup>/<sub>2</sub>



and support.

measurements

H1:20 + H2:7 cm 1+2

H1:30 + H2:19 cm



# Play box

Only for Rabbit. Play box with tray for different activities and play

size measurements

one H:9 x W:31 x D:10 cm



#### Extension piece Enlarges the centre bar and

allows even more precise

size measurements

one H:8 cm



#### Tray Only for Toucan. Comes as standard with activity bowl and cover plate

size measurements one D:61 x W:52 cm



#### Armrests for tray Only for tray for Toucan. The armrests support the forearms

size measurements

one W:25 x D:351/2 cm



# SNUG SERT

H:41/4"xW:6"xD:3"

size measurements