



FIT RIGHT GUIDE

See center insert page L

Dynamic Stander**WDS-150 Kaye Dynamic Stander – Early Intervention/Primary**

Age ranges:	18 mos. - 10 years
Pelvic block height (A):	20" - 29½"
Pelvic width (G):	5" - 11½"
Pelvic depth:	4½" - 8½"
Overall base dimensions:	24"W x 28"L
Tray dimensions:	26"W x 20"D
Casters:	4" swivel locking
Weight limit:	120 lbs.
Height range (E):	2'8" - 3'11"

WDS-150AB Abductor

Pad dimensions:	4"W x 7"H
Width adjustment:	4" - 6"

WDS-150LP Lateral Pads (pair)

Width adjustment:	5¼" - 12¼"
-------------------	------------

WDS-250 Kaye Dynamic Stander – Intermediate/Young Adult

Age ranges:	7 years - adult
Pelvic block height (A):	29½" - 45"
Pelvic width (G):	9" - 14½"
Pelvic depth:	5½" - 12½"
Overall base dimensions:	33"W x 33"L
Tray dimensions:	26"W x 22½"D
Casters:	4" swivel locking
Weight limit:	220 lbs.
Height range (E):	3' 11" - 6' 2"

WDS-250AB Abductor

Pad dimensions:	4"W x 7"H
Width adjustment:	4" - 6"

WDS-250LP Lateral Pads (pair)

Width adjustment:	8" - 13½"
-------------------	-----------

Marvel Vertical Stander**MV-100 Early Intervention/Preschool Marvel Vertical Stander**

User height (E):	25" - 41"
Tray height (C2):	18" - 29"
Hip width (G):	15"
Weight capacity:	75 lbs.
Base size:	22" x 23"

MV-200 Primary Marvel Vertical Stander

User height (E):	33" - 57"
Tray height (C2):	25½" - 43½"
Hip width (G):	15"
Weight capacity:	125 lbs.
Base size:	22" x 23"

NEW!

Kaye Dynamic Stander

The **Kaye Dynamic Stander** has all the features of a vertical stander, but in addition, the Dynamic Stander allows the client to shift weight within his base of support while standing. This permits the client to practice balance reactions that occur as individuals actively reach, bend, and move in standing. The Dynamic Stander helps develop steady-state balance while the client is safely aligned.

The Dynamic Standers, available in two sizes, are designed to accommodate children from 18 months through adulthood (up to 5' 10") who need support at the hips and pelvis to maintain standing alignment. Adjustable, lateral pelvic cushions can be used firmly against the sides of the pelvis to assure symmetry, while lower knee cushions can assist knee extension and help maintain hip-knee-foot alignment. Both standers come standard with a tray and swivel, locking casters. Optional accessories include an Abductor and Lateral Pads to help maintain hip-knee-foot alignment.

(Please note: availability for WDS-150 is late summer, 2007)



WDS-250

Marvel Vertical Stander

Simple and effective, the **Marvel Vertical Stander** provides vertical alignment and upright support. It is easy to assemble, adjust, and reposition.

The Marvel offers four points of support at the foot, knee, pelvic, and trunk level. It is available in two sizes. It comes standard with the following: foot restraint and heel cups with lateral and anterior/posterior adjustments to provide foot and ankle alignment; lateral, depth, and height adjustable knee blocks; a posterior pelvic stabilizer with a swing-open feature; height adjustable trunk support; removable, height and angle adjustable tray (18" x 23"); and base with 4 swivel wheels with locking casters.



MV-100

MV-200