

TherAdapt® Products Inc.

922 N. Washington Avenue, Ludington, MI 49431

Phone: 231-480-4008 • Fax: 231-843-3905

www.theradapt.com • customerservice@theradapt.com

MEASUREMENTS FOR THERADAPT SCHOOL CHAIRS

Model #	Seat depth (B1)	Seat width (G)*	Seat height (A) to floor	Seat height (A) to FPA-H	Seat height (A) to FPA-L	Back height (D)	Overall width	Armrest height to Seat (C1)	Armrest height to Floor (C2)
SC-100	9"-12½"	6½"-11½"	10"-14 ½"	3"-12 ½"	6"-12 ½"	12 ½"	17"	3 ½"-7 ½"	13"-22"
SCW-100	9"-12½"	8½"-12½"	10"-14 ½"	3"-12 ½"	6"-12 ½"	12 ½"	20"	3 ½"-7 ½"	13"-22"
SCX-100	9"-12½"	12"-17"	10"-14 ½"	4½"-14	7½"-14"	12 ½"	23"	3 ½"-7 ½"	13"-22"
SC-200	11½"-15"	8½"-13½"	13"-17 ½"	6"-15"	9"-15"	14"	20"	3 ½" – 11"	16 ½"-28 ½"
SCW-200	11½"-15"	12"-17"	13"-17 ½"	6"-15"	9"-15"	14"	23"	3 ½" – 11"	16 ½"-28 ½"
SCX-200	11½"-15"	15"-20"	13"-17 ½"	7½"-16½"	10½"-16½"	14"	26"	3 ½" – 11"	16 ½"-28 ½"
SC-300	13½"-16½"	12"-17"	14"-18 ½"	7"-16 ½"	10"-16 ½"	16"	23"	4 ½"-10 ½"	18 ½"-29"
SCW-300	13½"-16½"	15"-20"	14"-18 ½"	7"-16 ½"	10"-16 ½"	16"	26"	4 ½"-10 ½"	18 ½"-29"
SCX-300	13½"-16½"	18"-23"	14"- 18 ½"	7"-16 ½"	10"-16 ½"	16"	29"	4 ½"-10 ½"	18 ½"-29"

* Narrowed seat width achieved by using Side Pads.

MEASUREMENTS IN INCHES:

A _____ Lower leg length (measured from the back of the knee to the floor)

B1 _____ Upper leg length (measured from the back of the knee to the back of the seat as the child leans forward slightly)

C1 _____ Elbow height (measured from the bent elbow to the seat)

C2 _____ Elbow height (measured from the bent elbow to the floor)

D _____ Shoulder height (measured from the shoulder to the seat)

G _____ Hip width (measured at the widest location)