MOBILITY AIDS WWW.THERADAPT.COM

FIT RIGHT GUIDE See center insert page L

Posture Control Walkers

Height Adj. (B): Mid-buttock to floor **DESCRIPTION: 2 WHEELS**

14¹/₂" - 20" 16¹/₂" - 22" 19" - 25" 23" - 30¹/₂" 28" - 36" 35¹/₂" - 41¹/₂" W-1B W-2B* W-3B*

W-4B* W-5C* **DESCRIPTION: 4 WHEELS**

14¹/₂" - 18" 16¹/₂" - 22" W-1/2BR W-1BR W-2BR* 19" - 25" 23" - 30¹/₂" 28" - 36" 35¹/₂" - 41¹/₂" W-3BR* W-4BR* W-5CR*

DESCRIPTION: FRONT SWIVEL WHEELS

14¹/₂" - 18" 16¹/₂" - 22" 19" - 25" 23" - 30¹/₂" W-1/2BS W-1BS W-2BS* W-3BS* W-4BS* 28" - 36" 351/2" - 411/2" W-5CS*



The ORIGINAL reverse walker is available in three styles and five sizes. For all walkers, the base increases as the height of the walker is adjusted. They are available with three different wheel configurations: 2 Wheels with wheels in the front and walker tips in the back; 4 Wheels with ratcheted rear wheels

> that prevent the walker from rolling backwards during ambulation, and Front Swivel Wheels with the same ratcheted rear wheels but front swivel wheels that lock for straight-ahead tracking and unlock for easier turns. All

walkers fold compactly. In most cases Conversion Kits are available to convert the style of the walkers.



These walkers add the convenience of resting in a seated position whenever the user is tired and lacks the endurance for extended standing and walking, or wishes to sit in order to use both hands for play, self-care, or daily skills.

The seat folds up into the back space of the walker, so that it is out of the way for walking, and conveniently folds down for resting. The seat is also height adjustable, to better position the user with feet resting on the floor. The standard ratchet mechanism on the rear wheels prevent the walker from rolling backwards during ambulation. The walkers fold with the seat in place. They are available in all wheel combinations and can be used with any forearm support needed.



PostureRest Walkers

Height Adj. (B): Mid-buttock to floor

DESCRIPTION: 2 WHEELS WITH SEAT W-1/2BH

14¹/₂" - 20" 16¹/₂" - 22" W-1BH 20¹/₂" - 25¹/₂" 25¹/₂" - 31¹/₂" 30¹/₂" - 36" 35¹/₂" - 41¹/₂" W-2H* W-3H* W-4H* W-5CH* **DESCRIPTION: 4 WHEELS WITH SEAT** W-1/2BHR 14¹/₂" - 18" 16¹/₂" - 22"

W-1BHR 201/2" - 251/2 W-2HR* 25¹/₂" - 31¹/₂" 30¹/₂" - 36" 35¹/₂" - 41¹/₂" W-3HR* W-4HR* W-5CHR*

DESCRIPTION: FRONT SWIVEL WHEELS W/SEAT

W-1/2BHS 14¹/₂" - 18" 16¹/₂" - 22" W-1BHS 20¹/₂" - 25¹/₂" 25¹/₂" - 31¹/₂" 30¹/₂" - 36" 35¹/₂" - 41¹/₂" W-2HS* W-3HS* W-4HS* W-5CHS*

*Legs can switch to make a forward facing walker

Variable Resistance Rear Wheels



The Variable Resistance One-Way Ratchet Rear Wheels provide resistance to forward motion of the Kaye Posture Control and PostureRest walkers. These wheel kits can be retrofitted to all Kaye B, C, and H frame walkers.

Variable Resistance Rear Wheels

DESCRIPTION: REAR LEGS W/WHEELS

Model No. Compatible Walkers Model No. Compatible Walkers W2H, R, S (pair) VR1/2B W1/2B, H, R, S (pair) VR2H W3H, R, S (pair) W4H, R, S (pair) W5H, R, S (pair) W1B, R, S (pair) W2B, R, S (pair) W3B, R, S (pair) VR1B VR3H VR2B VR4H VR3R VR5H W4B, R, S (pair) VR4B VR5C W5C, H, R, S (pair)