





Booster Seat

Soft Touch® Booster Seat




child restraint


press center to unlock



insert tabs to lock



attachment strap




⚠ WARNING: SECURE ATTACHMENT STRAP TIGHTLY AROUND SEAT OF ADULT CHAIR. AUDIBLE “CLICK” ENSURES SAFE ATTACHMENT.

⚠ WARNING: ALWAYS SECURE CHILD IN THE RESTRAINT UNTIL CHILD IS ABLE TO GET IN AND OUT OF THE BOOSTER SEAT WITHOUT HELP (ABOUT 2.5 YEARS OLD). ALWAYS USE ATTACHMENT STRAP PROVIDED TO SECURE BOOSTER SEAT TO AN ADULT CHAIR. NEVER USE AS A BATH SEAT. THIS PRODUCT IS ONLY INTENDED FOR USE IN A CHAIR AS A BOOSTER. DO NOT USE ON ANY OTHER ELEVATED SURFACE.

⚠ WARNING: PREVENT SERIOUS INJURY OR DEATH. DO NOT USE IN MOTOR VEHICLES.

Specifications



Limited Warranty

Bergeron by Design warrants to the original retail purchaser of the Special Tomato® product, that if any part or component proves defective in material or workmanship within five years of the purchase date, the defective part will be repaired or replaced (at Bergeron by Design's discretion) free of charge. Warranty service may be performed by an authorized service center or (at Bergeron by Design's discretion) the factory.

This warranty does not cover normal wear and tear, or damage caused by accident or misuse.

To exercise this limited warranty, the user should first obtain a Return Authorization Number from the dealer the Special Tomato® product was purchased from. The product must be delivered charges pre-paid (UPS recommended) to the factory or to an authorized service center, together with a copy of the original invoice, the Return Authorization Number and a written description of the problem.

Product Care

Special Tomato® Products are impermeable to fluids. However, DO NOT SUBMERGE while cleaning. Use a mild detergent to clean periodically.

Customer Service

Contact your dealer directly. Or...

If your Special Tomato® product was purchased from the Adaptivemall.com® please use the following:

Adaptivemall.com® Customer Service Support

Toll Free: 1.800.371.2778

International: 1.315.429.7112

Fax: 1.315.429.8862

Email: info@adaptivemall.com

Notice

No part of this document may be photocopied, reproduced, transmitted, transcribed, stored in a retrieval system or translated to another language or computer language, in any form or by any means, electronic, mechanical, magnetic, optical, chemical, manual, or otherwise without the prior written consent of Bergeron by Design.

Use only Special Tomato® accessories and parts on Special Tomato® products. Special Tomato® parts are not interchangeable with other manufactures' products. Replace any worn parts immediately for safety. © Bergeron By Design.

! WARNING: Prevent serious injury or death from falls or sliding out

- Always secure child in the restraint
- Never leave child unattended
- Always keep child in view
- Always use attachment strap provided to secure to an adult chair



Bergeron By Design
15 South Second St.
Dolgeville, NY 13329
1.866.529.8407

