## **MODULAR CHAIRS**

## WWW.THERADAPT.COM


## Feeding Chair HEIGHT CHAIR


These two distinctly different chairs, the **Height Chair** and the **2-in-1 Booster Seat**, have both been designed with interaction, socialization, and age-appropriate mealtime

experiences in mind. They also provide a great place to play!

Appropriate postural support and positioning of the physically challenged child during mealtime is a key element in the ability to chew and swallow. The adjustability designed into the Height Chair and 2-in-1 Booster Seat allows for







## OPTIONS

- A. A Heel Stop, complete with dual ankle straps and curved side release buckles, is available to help keep the feet secure on the footplate.
- **B.** The **Height Chair** shown with the Abduction Block, Adjustable Armrest, Butterfly Kit, Foot Sandals, Headrest, Seat and Back Pad Set, Support Blocks (4), and Head Support Wedges. See pages 20-21.

FIT RIGHT GUIDE See center insert page L
Seat width G Back height P
Height ChairFC-100 Feeding Chair – Height ChairSeat depth (B1): $6" - 14"$ Seat width (G): $11^{1}/2"$ ( $6"$ to $11^{1}/2"$ )*Seat height to footrest (A): $4^{1}/2"$ to $15"$ Seat height to floor: $21" - 22 1/4"$ Side height to floor: $26"$ w/armrests: $27^{1}/2" - 30"$ Back height (D): $15"$ FC-100AB Abduction BlockFC-100AR Armrests (pair)Armrest height (C1): $6" - 9^{1}/2"$ Armrest height to floor: $27^{1}/2" \times 30"$
FC-100BK Butterfly Kit FC-100FS Foot Sandals (pair) Inside dimensions: 3"W x 5"L Outside dimensions: 41/4"W x 51/2"L FC-100HR Headrest Back height (E): 21" - 26" FC-100HS Heel Stop FC-100PS Seat and Back Pad Set FC-100PSW Seat (full width) and Back Pad Set FC-100SB Support Blocks* 5"H x 4"W (pair) FC-100SK Support Kit (includes FS, PS, and 2 pair SB)
FC-100SW Head Support Wedges 6"H x 6"W (pair)   FC-100T Wood Tray** Tray dimensions: 20"W x 18"D   Tray cutout: 14"W x 4"D   FC-100TC Clear Polycarbonate Tray**   Tray dimensions: 20"W x 18"D   Tray dimensions: 20"W x 18"D   Tray dimensions: 20"W x 18"D   Tray cutout: 13"/2"W x 4"D   *Note: Seat width reduction is achieved by use of Support Blocks.

\*Note: Seat width reduction is achieved by use of Support Blocks. \*\*Note: Armrests must be ordered to attach the Tray to the chair.