

## Corner Back Insert & Low Back Insert

The TherAdapt® **Corner Back Insert** and **Low Back Insert** interface with the three sizes of Adjustable Positioning Chairs and Bolster Chairs, allowing for seat depth reduction as needed. The inserts were designed to provide lateral and posterior trunk support for the user who does not require head support.

The **Corner Back Insert** provides thoracic support for the mild or moderately involved user who typically hyperextends at the trunk and retracts the scapulae. With the pelvis properly positioned in the appropriate chair, the insert encourages scapular protraction and normal spinal alignment. It comes standard with a mounting yoke, 3 Back Liner Pads, and a Butterfly Kit as shown.

The **Low Back Insert** provides thoracic support for mild or moderately involved users who typically flex forward or laterally at the trunk. With the pelvis properly positioned in the appropriate chair, the insert encourages thoracic extension and normal spinal alignment. It comes standard with a mounting yoke, 3 Back Liner Pads, a pair of thoracic pads, and an anterior chest strap as shown.




*Erin, age 5, is seated in the APC-200 supported by the LBI-200.*



*Joey, age 2 1/2, is seated in the Adjustable Positioning Chair (APC-200), with his hands in midline, supported by the Corner Back Insert (CBI-200).*

### FIT RIGHT GUIDE

See center insert page L



#### Corner Back Insert

##### CBI-100 – Early Intervention


Height (D): 11" - 15"  
Trunk width (F): 9"

##### CBI-200 – Preschool/Primary

Height (D): 12 1/2" - 18"  
Trunk width (F): 11"

##### CBI-300 – Primary/Intermediate

Height (D): 15" - 22"  
Trunk width (F): 14"



#### Low Back Insert

##### LBI-100 – Early Intervention

Height (D): 13 1/2" - 21"  
Trunk width (F): 5" - 9"

##### LBI-200 – Preschool/Primary

Height (D): 15 1/2" - 24"  
Trunk width (F): 7" - 11"

##### LBI-300 – Primary/Intermediate

Height (D): 17 1/2" - 28"  
Trunk width (F): 9" - 13"