MOBILITY AIDS WWW.THERADAPT.COM

KAYE PRODUCTS Anterior Chest Support Walkers

The Kaye Anterior Chest Support Walkers with or without forearm supports, have been specifically designed for users with severe spasticity, athetosis, or ataxia. They allow individuals who can support weight on their legs, but

lack the upper extremity strength to use a traditional anterior walker, the opportunity to stand, take steps, and move!

The front wheels can swivel or lock in the forward position while the rear wheels have an external ratchet mechanism to prevent the walker from rolling backwards. The anterior chest pad adjusts in height, width, and angle for individual customization. On the walker without forearm support the handgrip position adjusts front to back for comfortable, secure upper extremity positioning and support. Children who need additional support at the arms will find the adjustable forearm supports with multipositional handgrips ideal. The walkers fold flat for storage or transport.



ANTERIOR CHEST SUPPORT WALKER WITHOUT FOREARM SUPPORT

Suspension Accessory

The Kaye Suspension Accessory attaches to Kaye Posture Control B-series walkers and allows these walkers to function as a partial weight-bearing system. The Suspension Accessory, with the harness, holds the client securely and allows the client to practice over ground walking with partial or full weight-bearing as needed.

The walker folds with the Suspension Accessory in place for easy transportation or storage. The Suspension Accessory can be easily removed from the walker for more compact folding or when it is appropriate for the client to use the walker independent of the accessory.

Harnesses (sold separately) are available in three sizes. The small and medium slim harnesses are "diaper style" and fit between the legs of the child. The medium and large harnesses attach with thigh cuffs and circumferential support around the trunk and pelvis. Compression straps are included to control alignment or postural problems.



HARNESS WITH SUSPENSION **ACCESSORY**





Anterior Chest Support Walker WITHOUT FOREARM SUPPORT

Model No. Height Adj. (D)* Y-1S 22" - 36"* Y-25 30" - 47"* Y-35

*Note: This measurement is taken from axilla (arm pit) to the floor.

WITH FOREARM SUPPORT

Model No. Height Adj. (C1)** 19" - 26"** Y-1FS 211/2" - 30"** Y-2FS Y-3FS 28" - 39"** 37" - 45"**

**Note: This measurement is taken from the bent elbow to the floor.



SUSPENSION ACCESSORIES AND WALKERS

Suspension Accessory

Compatible Walker	Weight Limit
All W1/2B and W1B Walkers	60 lbs.
All W2B Walkers	85 lbs.
All W3B Walkers	130 lbs.
All W4B Walkers	180 lbs.
	All W1/2B and W1B Walkers All W2B Walkers All W3B Walkers

Harnesses for Suspension Accessory

Model	Size	Max. Pelvic Circum.	Weight Limit
WHS1	Small	16" - 25"	60 lbs.
WHS1.5	Medium Slim	18" - 28"	60 lbs.
WHS2	Medium	24" - 35"	200 lbs.
WHS3	Large	34" - 45"	200 lbs.