MOBILITY AIDS WWW.THERADAPT.COM



Miniwalk®

MIN-100 Miniwalk® - Early Intervention/Preschool

11/2 - 4 years Age ranges: 101/2" - 19" Inseam (B): 20" - 29" Trunk height (D): Overall dimensions: 251/2'W x 29"L 33 lbs. Weight: 90 lbs. Max user weight: Color. Turquoise

MIN-200 Miniwalk® - Preschool/Primary

4 - 7 years Age ranges: 14" - 25" Inseam (B): 27" - 35" Trunk height (D): Overall dimensions: 251/2"W x 29"L 33 lbs. Weight: Max user weight: 90 lbs. Color: Turquoise

MIN-AT Anti-tippers (pair) MIN-DB Drag Brakes (pair) MIN-LP Lateral Hip Pads (pair) MIN-LS Leg Separator Panel MIN-NR No Roll-Back Brakes (pair)

Meywalk®

MEY-300 Meywalk® - Intermediate

7 - 14 years Age ranges: 21" - 29" Inseam (B): Trunk height (D): 33" - 40" Overall dimensions: 271/2"W x 371/2"L Weight: 55 lbs. Max user weight: 220 lbs.

Blue Color: MEY-400 Meywalk® - Adult

14 years - Adult Age ranges: Inseam (B): 27" - 37" 41" - 51" Trunk height (D): Overall dimensions: 271/2"W x 371/2"L 57 lbs. Weight: Max user weight: 220 lbs.

Color: Blue MEY-AT Anti-tippers (pair) MEY-DB Drag Brakes (pair)

MEY-HL Height Lowering Kit (3") MEY-LP Lateral Hip Pads (pair)

MEY-LS Leg Separator Panel

MEY-NR No Roll-Back Brakes (pair)

ADULT MEYWALK



Meywalk® & Miniwalk® Gait Trainers

Miniwalk® from Meyland-Smith in Denmark is the newest model of walking supports for small children from approximately 1½ to 7 years of age. The height of the trunk support and seat are adjusted quickly with spanner grips. The seat and trunk support has a dynamic spring suspension to encourage a more natural gait pattern. The body posture of the child can be set in any position from vertical to 15° forward tilt, making it possible to find the body inclination that gives the child the best control. The Miniwalk® is available in two sizes, and both come standard with 6" front swivel wheels and 10" rear wheels.

The Intermediate and Adult Meywalk® Gait Trainer is the improved version of the famous MEYWALK® walking support. The height of the trunk support and seat are now adjusted quickly with spanner grips. The seat and trunk support has a dynamic spring suspension to encourage a more natural gait pattern. For transfers, the seat unit can be lowered by using a swing bar, making it easier for the user to get in and out of the walking support. Once the transfer is complete, the swing bar is pushed down and the user is lifted into the standing position. This lever lift ensures easier transfers for the caregiver and more comfort for the user. The body posture of the user can be set in any position from vertical to 15° forward tilt, making it possible to find the body inclination that gives the user the best control. Both sizes come standard with

8" front swivel wheels and 16" rear wheels.

The Miniwalk® and Meywalk® are made with high quality materials that guarantee long life, maximum reliability and a minimum of maintenance, making them suitable for indoor as well as outdoor use.

Both are easily adjustable, have all steel construction, a durable frame finish, and soft no-flat tires and casters. Both models come standard with seat, trunk support, rear stop, handle bar, fender wheels, and parking brakes. Available accessories include Anti-tippers,

> Pads, a Leg Separator Panel, and No-Roll-Back Brakes. In addition, a Height Lowering Kit is available for the Meywalk® that lowers the seat height 3".

Drag Brakes, Lateral Hip

PRESCHOOL/PRIMARY MINIWALK®